

USB® FidelityTaq™ DNA Polymerase

- Increases amplification fidelity approximately 3-6 times over Taq DNA Polymerase
- Allows for amplification of long PCR products
- Generates products whose ends are compatible with either blunt-end or TA cloning procedures

USB FidelityTaq DNA Polymerase

FidelityTaq DNA Polymerase combines high quality USB recombinant Taq DNA Polymerase with a highfidelity, proofreading polymerase. This enzyme blend has the 5'→3' exonuclease activity of Taq DNA polymerase as well as the 3'→5' exonuclease activity of the proofreading polymerase. FidelityTaq DNA Polymerase increases amplification fidelity approximately 3-6 times over Taq DNA Polymerase alone and allows for amplification of longer product sizes⁽¹⁻⁴⁾.

FidelityTaq DNA Polymerase generates PCR products whose ends are compatible with either blunt-end or TA cloning procedures⁽⁵⁾. Although a range of PCR product sizes can be amplified (Fig. 1), optimization may be necessary for targets longer than 10 kb from complex templates. For amplification of longer targets with less optimization and greater yields, use FidelityTaq PCR Master Mix or FidelityTaq PCR Master Mix Plus (PN 71182/71183).

Storage buffer:

20 mM Tris-HCl, pH 8.5, 1 mM DTT, 0.1 mM EDTA, 100 mM KCl, 50% glycerol, stabilizers.

Concentration:

5 units/μl

Functional tests:

Tested by amplifying a 20.7 kb PCR product from lambda DNA.

Functionally Tested 10X PCR Reaction Buffer (1 ml included, PN 71165):

100 mM Tris-HCl (pH 8.6), 500 mM KCl, 15 mM MgCl₂

Functionally Tested MgCl₂ (1 ml included, PN 71167):

25 mM solution

Purity:

Free from detectable non-specific nucleases.


FidelityTaq DNA Polymerase

Product code	Pack size	List price
71180	50 units	\$37.00
	250 units	\$161.00
	1,000 units	\$550.00
	5 × 250 units	\$821.00
	5,000 units	\$2,897.00

References:

- Barnes, W. M. (1994) *Proc. Natl. Acad. Sci. USA* **91**, 2216-2220.
- Cheng, S., Fockler, C., Barnes, W. M. and Higuchi, R. (1994) *Proc. Natl. Acad. Sci. USA* **91**, 5695-5699.
- Barnes, W. M. (1992) *Gene* **112**, 29-35.
- Cline, J., Braman, J. C. and Hogrefe, H. H. (1996) *Nucleic Acids Res.* **24**, 3546-3551.
- Magnuson, V. L., Ally, D. S., Nylund, S. J., Karanjawala, Z. E., Rayman, J. B., Knapp, J. I., Lowe, A. L., Ghosh, S. and Collins, F. S. (1996) *BioTechniques* **4**, 700-709.
- Innis, M. A., Myambo, K. B., Gelfand, D. H. and Brow, M. A. (1988) *Proc. Natl. Acad. Sci.* **85**, 9436-9440.
- Lawyer, F. C., Stoffel, S., Saiki, R. K., Myambo, K., Drummond, R. and Gelfand, D. H. (1989) *J. Biol. Chem.* **264**, 6427-6437.

Fig. 1. Amplification of diverse targets



Target (source): Lane 1: Numb, 455 bp (10 ng human genomic DNA); Lane 2: NRAGE, 1.55 kb (10 ng human genomic DNA); Lane 3: Numb, 4.6 kb (10 ng human genomic DNA); Lane 4: Lambda target, 20.7 kb (10 ng lambda DNA). A total of 25 cycles was used for lambda target and 35 cycles for human genomic DNA targets. A final concentration of 1.5 mM MgCl₂ was used in all reactions.

Related products:

Amplify

HotStart-IT® Taq DNA Polymerase [71195]

Includes 10X PCR Buffer and 25 mM MgCl₂
50 units
250 units
1,000 units
5 × 250 units
5,000 units

HotStart-IT Taq Master Mix (2X) [71196]

25 reactions
100 reactions
500 reactions

RubyTaq™ DNA Polymerase [71190]

Includes 10X PCR Buffer and 25 mM MgCl₂
50 units
250 units
1,000 units
5,000 units

RubyTaq PCR Master Mix (2X) [71191]

100 reactions (50 µl/rctn)
500 reactions (50 µl/rctn)

USB® Taq DNA Polymerase [74165]

Includes 10X PCR Reaction Buffer
50 units
400 units
2,000 units
4,000 units (2 × 2,000 units)
20,000 units

Taq PCR Kit [71161]

100 reactions (100 µl/rctn)

Taq PCR Master Mix, 2X [71162]

100 reactions (50 µl/rctn)
4 × 625 µl per pack

Taq Master Mix Plus [71163]

Contains mix, 25 mM MgCl₂ and PCR-Qualified H₂O.
100 reactions (50 µl/rctn)

PCR Nucleotide Mix, 10 mM [77212]

500 µl
2 × 500 µl

PCR Nucleotide Mix, 25 mM [77119]

500 µl

dNTP Set of Four (dA, dC, dG, dT) 100 mM solutions [77100]

4 × 25 µmol (250 µl)
4 dNTPs per pack

Betaine, 5 M Solution [77507]

1.5 ml
5 × 1.5 ml
10 ml

Magnesium Chloride, 1 M Solution [78641]

10 × 1 ml
100 ml

Mineral Oil [71600]

10 ml
25 ml
1 L

Purify

ExoSAP-IT Express PCR Product Cleanup [75001]

20 reactions
100 reactions
500 reactions
2,000 reactions
5,000 reactions
480 reactions × 8-tube strip

HT ExoSAP-IT Fast High-Throughput PCR Product Cleanup [78595]

20 reactions
480 reactions × 8-tube strip
5,760 reactions (12 × 8-tube strips in a tray)
23,040 reactions (48 × 8-tube strips in a tray)
1,000 reactions (2 ml)
5,000 reactions (10 ml)
1,920 reactions (96-well plate, 40 µl per well)
11,520 reactions (6 × 96-well plate, 40 µl per well)

Biochemicals

Agarose - Low Melt, Separation, ≤1000 bp [32829]

Genetic Performance Certified™
25 gm
100 gm

Agarose - Separation ≥ 500 bp [75817]

Genetic Performance Certified
25 gm
100 gm
250 gm
500 gm

Agarose - LE [32802]

25 gm
100 gm
250 gm
500 gm
1 kg

TBE Buffer, 5X Solution [75891]

1 L
5 L

TE Buffer, 1X Solution [75893]

10 × 1 ml
100 ml
500 ml

Water, Nuclease-Free [71786]

10 × 1 ml
100 ml
500 ml
1 L
5 L

For Research Use Only. Not for use in diagnostic procedures.

Affymetrix, USB, ExoSAP-IT, and HotStart-IT are registered trademarks of Affymetrix, Inc. FidelityTaq, Genetic Performance Certified and RubyTaq are trademarks of Affymetrix, Inc. Please see affymetrix.com/trademarks for a complete list of Affymetrix trademarks. ExoSAP-IT PCR Product Cleanup product line—These products are covered by patents owned by Affymetrix, including US Patent Nos. 6,379,940 and 6,387,634 and corresponding patents in Canada and Japan. Taq DNA Polymerase—Purchase of this product conveys to the purchaser the limited, non-transferable right under certain patents to use only the purchased amount of product for the purchaser's own internal research. No right under any other patent claims or for any other use is conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting Out-licensing at Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad, California 92008; outlicensing@lifetech.com.

Affymetrix, Inc.

USB® Products

26111 Miles Road
Cleveland, Ohio 44128
Tel: 800-321-9322 | 216-765-5000
Fax: 800-535-0898 | 216-464-5075
USBcustomerserv@affymetrix.com
usb.affymetrix.com

Affymetrix UK Ltd.

USB® Products

Voyager, Mercury Park,
Wycombe Lane, Wooburn Green,
High Wycombe HP10 0HH, UK
Tel: +44 (0)1628 55 2600
Fax: +44 (0)1628 55 2675
USBcustomerserveurope@affymetrix.com

Affymetrix Pte Ltd.

USB® Products

7 Gul Circle, #2M-01
Keppel Logistics Building
Singapore 629563
Tel: +65 63957310
Fax: +65 63957300
USBcustomerserv@affymetrix.com