

PCR Tools

# USB® HotStart-IT® for real-time PCR

### Reliable performance

- Functionally tested to deliver a linear dynamic range of seven orders of magnitude
- USB HotStart-IT Taq DNA Polymerase only needs a 2 minute heat activation step
- Detect single copy targets
- Matches performance of popular master mixes such as ABI TaqMan® Universal PCR Master Mix

### Value and convenience

- Price per reaction is 25% less than TaqMan Universal PCR Master Mix
- Multiple formulations available for different detection chemistries

USB HotStart-IT master mixes for real-time PCR enable robust quantitative detection of nucleic acids. Designed to provide consistent results, HotStart-IT qPCR master mixes perform equally well as traditional real-time PCR reagents from other suppliers (Fig. 1) at a better price per reaction. And, unlike chemically modified Taq DNA Polymerase, USB HotStart-IT Taq DNA Polymerase only needs a 2 minute heat activation step. HotStart-IT qPCR master mixes are functionally tested to ensure linear performance over a broad dynamic range of seven orders of magnitude (Fig. 2).


USB HotStart-IT qPCR reagents utilize a novel primer binding protein to achieve a hot start. Primers are sequestered at lower temperatures which prevent the formation of non-specific products and primer dimers. The HotStart-IT binding protein enables efficient and specific amplification of the desired product, resulting in more sensitive real-time PCR. In addition, HotStart-IT master mixes only need a 2 minute heat activation step in contrast to chemically modified polymerases that require excessive heat activation steps over 10 minutes.

Formulations are available for either specific (e.g. TaqMan) or non-specific (e.g. SYBR®) detection chemistries. Mixes containing Uracil-DNA Glycosylase are also available if carry-over contamination prevention is desired. Additionally, HotStart-IT qPCR reagents are available with the passive reference dye, ROX™, already mixed in the master mix, or with both ROX and fluorescein included in separate tubes for maximum flexibility.


For expression analysis of RNA using qRT-PCR, HotStart-IT qRT-PCR master mixes are available in convenient one-step formats to perform both cDNA synthesis and real-time PCR in the same tube. If a two-step qRT-PCR workflow is desired, the USB First Strand cDNA Synthesis Kit for Real-Time PCR can be used to reliably generate cDNA prior to qPCR. This kit eliminates end bias that can be introduced by traditional reverse transcription methods utilizing either oligo(dT) or random hexamers.

Try USB HotStart-IT master mixes and maximize the value of your real-time PCR.

**Fig. 1. Performance comparison of HotStart-IT Probe qPCR Master Mix (PN 75766)**



**Fig. 2. Real-time PCR Amplification using HotStart-IT Probe qPCR Master Mix Kit (PN 75766)**



GAPDH Assay using HotStart-IT Probe qPCR Master Mix (PN 75766). Triplicate reactions were performed with a cloned region of the human GAPDH gene as template using an ABI 7500 Fast instrument. A TaqMan probe with FAM as the reporter fluorophore and BHQ-1<sup>®</sup> as the quencher was used to detect the 122 bp amplicon. ROX was used as a passive reference dye. The amplification process was linear over eight orders of magnitude (see inset) and a single copy of the target could be efficiently detected. The No Template Control (NTC) reaction generated no measurable fluorescence.

**HotStart-IT real-time PCR master mixes**

	100 x 50 µl reactions	500 x 50 µl reactions
HotStart-IT Probe qPCR Master Mix (2X), PN 75766		
HotStart-IT SYBR Green qPCR Master Mix (2X), PN 75762		
HotStart-IT Probe qPCR Master Mix with UDG (2X), PN 75764		
HotStart-IT SYBR Green qPCR Master Mix with UDG (2X), PN 75760		

**HotStart-IT one-step qRT-PCR master mixes**

	100 x 50 µl reactions	500 x 50 µl reactions
HotStart-IT Probe One-Step qRT-PCR Master Mix Kit, PN 75772		
HotStart-IT SYBR Green One-Step qRT-PCR Master Mix Kit, PN 75770		

**Reverse Transcription - cDNA Synthesis**

	50 reactions
First Strand cDNA Synthesis Kit for Real-Time PCR, PN 75780	

For research use only. Not for use in diagnostic procedures. Affymetrix, USB and HotStart-IT are registered trademarks of Affymetrix, Inc. Tested User Friendly is a trademark of Affymetrix, Inc. HotStart-IT Taq DNA Polymerase—Methods for using this product may be covered by US Patent No. 7,700,281. Taq DNA Polymerase—sold under licensing arrangements with Applied Biosystems. Purchase is accompanied by a limited license to use it in the Polymerase Chain Reaction (PCR) process in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Perkin-Elmer or as purchased, i.e., an authorized thermal cycler. TaqMan is a registered trademark of Roche Molecular Systems, Inc. ROX is a trademark of Applied Biosystems or its subsidiaries in the US and certain other countries. SYBR is a registered trademark of Molecular Probes, Inc. and is provided under an agreement with Molecular Probes, Inc. BHQ-1 is a registered trademark of Biosearch Technologies.

<p><b>Affymetrix, Inc.</b>  <b>USB® Products</b>                      26111 Miles Road                      Cleveland Ohio 44128                      Tel: 800-321-9322   216-765-5000                      Fax: 800-535-0898   216-464-5075                      USBcustomerserv@affymetrix.com  <a href="http://usb.affymetrix.com">usb.affymetrix.com</a></p>	<p><b>Affymetrix UK Ltd.</b>  <b>USB® Products</b>                      Voyager, Mercury Park,                      Wycombe Lane, Wooburn Green,                      High Wycombe HP10 0HH, United Kingdom                      Tel: +44 (0)1628 55 2600                      Fax: +44 (0)1628 55 2675                      USBcustomerserveurope@affymetrix.com</p>	<p><b>Affymetrix Pte Ltd</b>  <b>USB® Products</b>                      7 Gul Circle, #2M-01                      Keppel Logistics Building                      Singapore 629563                      Tel: +65 63957301                      Fax: +65 63957300                      USBcustomerserv@affymetrix.com</p>
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------