



Affymetrix® and Ambion® HT WT Expression Assay Consumables for Biomek® FX^P Target Prep Express

Intended Use

The Affymetrix® and Ambion® HT WT Expression Assay Consumables for Biomek® FX^P Target Prep Express Instrument contains plastic consumables needed to run target preparation on the Biomek® FX^P using the Ambion® WT Expression Kit combined with the Affymetrix® GeneChip® WT Terminal Labeling Kit. The accessory kit contains consumables sufficient for 4x24 or 4x96 rxn runs. Please refer to the *GeneChip® HT WT User Guide* (P/N 702923) for instructions on using the consumables.

Component	Number per Run	Number per Kit
96-Well Deep Well Microplate, polypropylene, 2 mL, non-sterile	2	8
Greiner U-bottom (round), polypropylene	2	8
BioRad Hard Shell PCR plate, blue shell/clear well	18	72
Greiner UV Star, Non-Sterile, clear, flat bottom	2	8
Costar polystyrene 96-well, U-bottom plate	4	16
24 well plate, 10 mL, square well, pyramid bottom	1	4
Pad, Microseal 'P'	1	2

Storage

The Affymetrix and Ambion HT WT Expression Assay Consumables for Biomek FX^P Target Prep Express should be stored at room temperature.

Ordering Information

P/N	Product Name	Description
901646	Affymetrix® and Ambion® HT WT Expression Assay Consumables for Biomek® FX ^P Target Prep Express	Labware sufficient for 4x24 or 4x96 rxn runs

Affymetrix® products can be purchased directly from Affymetrix in the United States, and many European and Asian countries. For all other territories, please view a list of our distribution partners, which can be located at: www.affymetrix.com/site/contact/index.affx.

Figure 1. 96-Well Deep Well Microplate




Figure 2. Greiner U-Bottom (Round), Polypropylene Plate




Figure 3. Bio-Rad Hardshell PCR Plate




Figure 4. Greiner UV Star Plate




Figure 5. Costar polystyrene 96-well, U-bottom plate





Figure 6. Plate, 24 well plate, 10 mL, square well, pyramid bottom



Figure 7. Auto-sealing Microplate Lid with Microseal P Pad



Contact Information

Affymetrix, Inc.

3420 Central Expressway
Santa Clara, CA 95051 USA
E-mail: support@affymetrix.com
Tel: 1-888-362-2447 (1-888-DNA-CHIP)
Fax: 1-408-731-5441

Affymetrix UK Ltd

Voyager, Mercury Park,
Wycombe Lane, Wooburn Green,
High Wycombe HP10 0HH
United Kingdom
E-mail: supporteurope@affymetrix.com
UK and Others Tel: +44 (0) 1628 552550
France Tel: 0800919505
Germany Tel: 01803001334
Fax: +44 (0) 1628 552585

Affymetrix Japan, K. K.

ORIX Hamamatsucho Bldg, 7F
1-24-8 Hamamatsucho, Minato-ku
Tokyo 105-0013, Japan
E-mail: supportjapan@affymetrix.com
Tel: +81-3-6430-4020
Fax: +81-3-6430-4021

Please visit our web site for international distributor contact information
www.affymetrix.com

Limited License

Subject to the Affymetrix terms and conditions that govern your use of Affymetrix products, Affymetrix grants you a non-exclusive, non-transferable, non-sublicensable license to use this Affymetrix product only in accordance with the manual and written instructions provided by Affymetrix. You understand and agree that except as expressly set forth in the Affymetrix terms and conditions, that no right or license to any patent or other intellectual property owned or licensable by Affymetrix is conveyed or implied by this Affymetrix product. In particular, no right or license is conveyed or implied to use this Affymetrix product in combination with a product not provided, licensed or specifically recommended by Affymetrix for such use.

Trademarks

Affymetrix®, Axiom™, Command Console®, DMET™, GeneAtlas™, GeneChip®, GeneChip-compatible™, GeneTitan®, Genotyping Console™, NetAffx®, and Powered by Affymetrix™ are trademarks or registered trademarks of Affymetrix Inc. All other trademarks are the property of their respective owners.

Copyright

© 2010 Affymetrix, Inc. All rights reserved.