

# Fetal Bovine Serum

From USB – Your Trusted Source

- Quality
- Purity
- Reliability



Fetal Bovine Serum (FBS) from USB is an excellent supplement for cell culture medium to properly grow and maintain cultured cells.

FBS from USB is provided as a sterile, cell culture-tested liquid derived from bovine blood that is collected and processed in the United States.

## Fetal Bovine Serum (FBS)

**Source:**  
Bovine Blood

## Product Specifications

**Form:**  
Clear straw to amber liquid

**Total Protein:**  
3.0 – 4.0 gm/dl

**Endotoxin (By LAL Gel Clot):**  
≤ 10 EU/ml

**Hemoglobin:**  
≤ 20 mg/dl

**pH:**  
6.7 – 8.0

**Mycoplasma (By Barile Method):**  
None detected

**Testing:**  
Tested and found negative for Bovine Viral Diarrhea (BVD), Infectious Bovine Rhinotracheitis (IBR) and Parainfluenza Type 3 Virus (PI3).

**Sterility:**  
Sterile

**Storage & Shipping:**  
Store at -20°C. Shipped on dry ice.

Product Code	Pack Size
15762	100 ml
15762	500 ml



©2004 USB Corporation. USB and logo design are registered trademarks of USB Corporation. The phrase 'Fueling Innovation' is a trademark of USB Corporation.

All goods and services sold are subject to the terms & conditions of sale from USB Corporation or the company which supplies them. A copy of these terms & conditions is available upon request.

**USB Corporation**  
26111 Miles Road; Cleveland, OH 44128  
800.321.9322 | [www.usbweb.com](http://www.usbweb.com)

 **usb**<sup>®</sup>  
Fueling Innovation™